


# USB2 - CMOS 1024 - 500

Complete medium speed, low noise cmos linescan camera system



## Key Features:

- NMOS-line scan camera system.
- 256 to 1024 Pixels low noise.
- 12 Bit ADC
- 4.000:1 Dynamic Range
- 2.2 ms to > 100 s exposure time
- USB 2.0 Interface, bus powered
- Drivers for Win 98 / XP / Labview.

## Overview:

The USB2-CMOS is an easy to use, complete ccd linescan camerasystem. It includes a low noise CMOS linescan camera with USB 2.0 interface. Additional components are not required.

The USB2-CMOS was designed for medium speed applications with the need for low noise precision light measurements.

## Applications:

- Spectroscopy.
- Portable applications

## Hardware:

The USB2-CMOS camera head includes the complete CMOS timing with signal conditioning (CDS), a precision 12 Bit ADC and an USB 2.0 interface.

The camera head is powered by the USB-bus. Additional power-supplies are not required. The USB2-CMOS provides optional start of scan output and an optional input for external triggering.

khs instruments

## Software

The USB2-CMOS linescan camera system is shipped with a software and drivers for Windows 98 and Windows XP.

The software includes a DLL to provide an interface to other software and the user software. Drivers for Labview are available upon request.

The DLL configures the USB2-CMOS by reading internal stored EEPROM data, so in most cases there is no need to configure the camera.




One single DLL supports all and up to 127 (different) khs-instruments cameras connected to the USB-bus.

The user software supports all cameras from khs-instruments. It provides following and more functions:

- X/Y Scaling
- X/Y Zoom
- Two Cursors
- Averaging
- Binning
- Subtract a Reference Scan

It provides functions to read, write and print stored (ASCII) files.

## Mechanical Dimension



## Specifications

### Detector Array:

Sensor: Hamamatsu S8378  
 Number of pixels: 256 to 1024  
 Pixel size: 25  $\mu\text{m}$  x 500  $\mu\text{m}$   
 Spectral range: 200 nm..1000 nm  
 Sensitivity nonuniformity: 3%  
 Saturation exposure Gain low: 0.57 lx sec

### System & Detector:

ADC resolution: 12 bits  
 Readout Noise Gain low: typ. 1 counts rms  
 Readout Noise Gain high: typ. 2 counts rms  
 Exposure Time: 2.2 ms to >5 s.  
 Linerate: 450 lps.

### System Requirements:

Operating system: Win 98 / XP.  
 Disk: 300 KB free.


### USB Interface:

USB: 2.0.  
 Required Current: 500 mA.

**Optional:**  
 Trigger: Input TTL.  
 Start of scan: Output TTL.

### Software:

Software includes: User software,  
 DLL interface,  
 Driver for Labview  
 upon request.



### User Software:

X scale edit: Enter start and stop.  
 Enter the values at two cursor positions.

Y scale edit: Enter start and stop.  
 Enter the values at two cursor positions.

X / Y unit edit: Enter units.

Averaging: Integration of several scans (up to 15).  
 Running mean of n consecutive scans.

Binning: Up to 64 pixels.

Display options: Display actual scan.  
 Load reference from actual scan and display scan minus reference.  
 Set reference to zero.

Data operations: Write to disk.  
 Write consecutive scans to disk.  
 Read from disk.  
 Print scan.